

ABX Minocal


LOT
CX 477
CAL

(Exp.) 2023-01-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	8.8			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.59	4.60			4.60					± 0.06
HB HGB	g/dl	14.0	14.0			14.0					± 0.2
	g/l	140	140			140					± 2
	mmol/l	8.69	8.69			8.69					± 0.12
HT HCT	%	39.9	39.8			40.3					± 1.0
	l/l	0.399	0.398			0.403					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	252	252			250					± 10
VMP MPV	µm ³ ; fl	8.9	8.6			8.0					± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300090308