

ABX Minocal

LOT		CX 473		CAL		(Exp.) 2022-09-05 (YYYY - MM - DD)									
		Rev 1		Minilysebio											
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE				
		ES60	CRP 200			60 8P 60 18P									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.1			9.2					± 0.2				
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.45			4.45					± 0.06				
HB HGB	g/dl	13.3	13.1			13.3					± 0.2				
	g/l	133	131			133					± 2				
	mmol/l	8.26	8.15			8.26					± 0.12				
HT HCT	%	37.5	38.0			38.5					± 1.0				
	l/l	0.375	0.380			0.385					± 0.010				
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	251			253					± 10				
VMP MPV	µm ³ ; fl	8.4	8.2			8.6					± 0.5				

ABX Minocal

LOT		CX 473		CAL		(Exp.) 2022-09-05 (YYYY - MM - DD)							
		Rev 1		LyseVet									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS								TOLERANCES TOLERANCE	
		Scil Vet ABC	Scil Vet ABC Plus ESV 60										
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.8									± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.42									± 0.06	
HB HGB	g/dl	13.0	13.0									± 0.2	
	g/l	130	130									± 2	
	mmol/l	8.07	8.07									± 0.12	
HT HCT	%	37.4	37.3									± 1.0	
	l/l	0.374	0.373									± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	240	250									± 10	
VMP MPV	μm ³ ;fl	N/A	8.7									± 0.5	

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208