

ABX Minocal

LOT

CX 474

CAL



(Exp.) 2022-10-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.1			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.55			4.45					± 0.06
HB HGB	g/dl	13.3	13.1			13.3					± 0.2
	g/l	133	131			133					± 2
	mmol/l	8.26	8.15			8.26					± 0.12
HT HCT	%	37.4	38.0			38.3					± 1.0
	l/l	0.374	0.380			0.383					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	248			248					± 10
VMP MPV	µm ³ ; fl	8.3	8.3			8.8					± 0.5

ABX Minocal

LOT

CX 474

CAL



(Exp.) 2022-10-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.45							± 0.06
HB HGB	g/dl	13.2	13.1							± 0.2
	g/l	132	131							± 2
	mmol/l	8.20	8.14							± 0.12
HT HCT	%	37.9	37.6							± 1.0
	l/l	0.379	0.376							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	237	240							± 10
VMP MPV	µm ³ ; fl	N/A	8.8							± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208