

ABX Minocal

LOT		CX 476		CAL		(Exp.) 2022-12-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.55			4.54					± 0.06
HB HGB	g/dl	13.4	13.0			13.4					± 0.2
	g/l	134	130			134					± 2
	mmol/l	8.32	8.09			8.32					± 0.12
HT HCT	%	37.8	37.8			38.6					± 1.0
	l/l	0.378	0.378			0.386					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	231	245			244					± 10
VMP MPV	µm ³ ; fl	9.3	8.9			9.8					± 0.5

ABX Minocal

LOT		CX 476		CAL		(Exp.) 2022-12-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	9.0								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.50								± 0.06
HB HGB	g/dl	13.0	13.0								± 0.2
	g/l	130	130								± 2
	mmol/l	8.07	8.07								± 0.12
HT HCT	%	37.1	37.6								± 1.0
	l/l	0.371	0.376								± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234	240								± 10
VMP MPV	µm ³ ; fl	N/A	9.6								± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208