

# ABX Minocal

**LOT**

**CX 473**

**CAL**



**(Exp.) 2022-09-05**  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.4</b>	<b>9.1</b>			<b>9.0</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.44</b>	<b>4.45</b>			<b>4.45</b>					± 0.06
HB HGB	g/dl	<b>13.2</b>	<b>13.4</b>			<b>13.3</b>					± 0.2
	g/l	<b>132</b>	<b>134</b>			<b>133</b>					± 2
	mmol/l	<b>8.20</b>	<b>8.32</b>			<b>8.26</b>					± 0.12
HT HCT	%	<b>37.5</b>	<b>38.0</b>			<b>37.8</b>					± 1.0
	l/l	<b>0.375</b>	<b>0.380</b>			<b>0.378</b>					± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>248</b>	<b>251</b>			<b>250</b>					± 10
VMP MPV	μm <sup>3</sup> ; fl	<b>8.2</b>	<b>8.2</b>			<b>7.4</b>					± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300090308