

ABX Minocal



(Exp.) 2023-01-05

(YYYY - MM - DD)

LOT		CX 477		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	8.8			9.1						± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.59	4.60			4.60						± 0.06
HB HGB	g/dl	14.0	14.0			14.0						± 0.2
	g/l	140	140			140						± 2
	mmol/l	8.69	8.69			8.69						± 0.12
HT HCT	%	39.9	39.8			40.3						± 1.0
	l/l	0.399	0.398			0.403						± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	252	252			250						± 10
VMP MPV	µm ³ ; fl	8.9	8.6			8.0						± 0.5

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.48