

ABX Minocal

LOT
CX 474
CAL

(Exp.) 2022-10-05

 Rev **1**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.1			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.55			4.45					± 0.06
HB HGB	g/dl	13.3	13.1			13.3					± 0.2
	g/l	133	131			133					± 2
	mmol/l	8.26	8.15			8.26					± 0.12
HT HCT	%	37.4	38.0			38.3					± 1.0
	l/l	0.374	0.380			0.383					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	248			248					± 10
VMP MPV	µm ³ ; fl	8.3	8.3			8.8					± 0.5

ABX Minocal

LOT		CX 474		CAL		(Exp.) 2022-10-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	9.0								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.45								± 0.06
HB HGB	g/dl	13.2	13.1								± 0.2
	g/l	132	131								± 2
	mmol/l	8.20	8.14								± 0.12
HT HCT	%	37.9	37.6								± 1.0
	l/l	0.379	0.376								± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	237	240								± 10
VMP MPV	µm ³ ; fl	N/A	8.8								± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208