

ABX Minocal

LOT		CX 475		CAL		(Exp.) 2022-11-05 (YYYY - MM - DD)					
Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2			9.4					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.59			4.49					± 0.06
HB HGB	g/dl	13.5	13.3			13.6					± 0.2
	g/l	135	133			136					± 2
	mmol/l	8.38	8.28			8.45					± 0.12
HT HCT	%	38.0	38.6			38.6					± 1.0
	l/l	0.380	0.386			0.386					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233	245			240					± 10

ABX Minocal

VMP MPV	μm^3 ; fl	8.2	8.3			8.6					± 0.5
---------	----------------------	-----	-----	--	--	-----	--	--	--	--	-----------

LOT	CX 475	CAL	 (Exp.) 2022-11-05 (YYYY - MM - DD)
	Rev 1		

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	$10^3/\text{mm}^3; 10^9/\text{l}$	9.1	9.1							± 0.2
GR RBC	$10^6/\text{mm}^3; 10^{12}/\text{l}$	4.51	4.53							± 0.06
HB HGB	g/dl	13.2	13.2							± 0.2
	g/l	132	132							± 2
	mmol/l	8.20	8.20							± 0.12
HT HCT	%	37.7	38.1							± 1.0
	l/l	0.377	0.381							± 0.010

48 FRONT / RECTO 1300105208

ABX Minocal



PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	218	240								± 10
VMP MPV	µm ³ ; fl	N/A	8.6								± 0.5

Ref: TEMP-0387 Rev. .