

ABX Minocal

LOT

CX 477

CAL



(Exp.) 2023-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7	8.8			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.60			4.60					± 0.06
HB HGB	g/dl	14.3	13.7			14.2					± 0.2
	g/l	143	137			142					± 2
	mmol/l	8.88	8.52			8.82					± 0.12
HT HCT	%	39.8	39.8			40.5					± 1.0
	l/l	0.398	0.398			0.405					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	240	252			254					± 10
VMP MPV	µm ³ ; fl	9.3	8.6			9.0					± 0.5

ABX Minocal

LOT

CX 477

CAL



(Exp.) 2023-01-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.60	4.58							± 0.06
HB HGB	g/dl	13.9	13.9							± 0.2
	g/l	139	139							± 2
	mmol/l	8.63	8.63							± 0.12
HT HCT	%	40.0	39.8							± 1.0
	l/l	0.400	0.398							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	224	250							± 10
VMP MPV	µm ³ ; fl	N/A	9.5							± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208