

ABX Minocal



LOT CX 473
Rev 1

CAL

(Exp.) 2022-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.1			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.45			4.45					± 0.06
HB HGB	g/dl	13.3	13.1			13.3					± 0.2
	g/l	133	131			133					± 2
	mmol/l	8.26	8.15			8.26					± 0.12
HT HCT	%	37.5	38.0			38.5					± 1.0
	l/l	0.375	0.380			0.385					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	251			253					± 10
VMP MPV	µm ³ ; fl	8.4	8.2			8.6					± 0.5

ABX Minocal



LOT CX 473

CAL

(Exp.) 2022-09-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	8.8	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.50	4.42							± 0.06
HB HGB	g/dl	13.0	13.0							± 0.2
	g/l	130	130							± 2
	mmol/l	8.07	8.07							± 0.12
HT HCT	%	37.4	37.3							± 1.0
	l/l	0.374	0.373							± 0.010
PLA PLT	10 ³ /mm ³ ; 10 ⁹ /l	240	250							± 10
VMP MPV	μm ³ ; fl	N/A	8.7							± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208