

ABX Minocal


LOT
CX 475
CAL

(Exp.) 2022-11-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	9.5	9.5	9.6	9.4	9.5	9.2	9.5	9.4	9.4	± 0.2
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.51	4.51	4.48	4.51	4.54	4.48	4.54	4.54	4.51	± 0.06
HB HGB	g/dl	13.3	13.3	13.3	13.4	13.4	13.3	13.4	13.4	13.4	± 0.2
	g/l	133	133	133	134	134	133	134	134	134	± 2
	mmol/l	8.26	8.26	8.26	8.32	8.32	8.26	8.32	8.32	8.32	± 0.12
HT HCT	%	37.9	37.7	36.7	36.8	36.8	36.3	36.8	36.8	36.8	± 1.0
	l/l	0.379	0.377	0.367	0.368	0.368	0.363	0.368	0.368	0.368	± 0.010
PLA PLT	10 ³ /mm ³ ; 10 ⁹ /l	242	249	240	242	245	237	245	244	242	± 10
VMP MPV	µm ³ ; fl	8.2	8.4	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

 1300090308
 Ref. TEMP-0387 Rev.48 BACK / VERSO