

ABX Minocal



(Exp.) 2022-10-05

(YYYY - MM - DD)

LOT		CX 474		CAL		(Exp.) 2022-10-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.1			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.55			4.48					± 0.06
HB HGB	g/dl	13.3	13.4			13.3					± 0.2
	g/l	133	134			133					± 2
	mmol/l	8.26	8.32			8.26					± 0.12
HT HCT	%	37.9	38.0			37.9					± 1.0
	l/l	0.379	0.380			0.379					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	247	248			242					± 10
VMP MPV	µm ³ ; fl	8.3	8.3			7.3					± 0.5

 1300090308
 Ref: TEMP-0387 Rev.48 FRONT / RECTO