

ABX Minocal



(Exp.) 2022-11-05

(YYYY - MM - DD)

LOT		CX 475		CAL		(Exp.) 2022-11-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.6	9.2			9.4					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.59			4.53					± 0.06
HB HGB	g/dl	13.4	13.6			13.4					± 0.2
	g/l	134	136			134					± 2
	mmol/l	8.32	8.45			8.32					± 0.12
HT HCT	%	38.0	38.6			38.3					± 1.0
	l/l	0.380	0.386			0.383					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	240	245			238					± 10
VMP MPV	µm ³ ; fl	8.3	8.3			7.5					± 0.5

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.48

ABX Minocal


 (Exp.) 2022-11-05
 (YYYY - MM - DD)

LOT		CX 475				CAL		(Exp.) 2022-11-05 (YYYY - MM - DD)			
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.5	9.5	9.6	9.4	9.5	9.2	9.5	9.4	9.4	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.51	4.48	4.51	4.54	4.48	4.54	4.54	4.51	± 0.06
HB HGB	g/dl	13.3	13.3	13.3	13.4	13.4	13.3	13.4	13.4	13.4	± 0.2
	g/l	133	133	133	134	134	133	134	134	134	± 2
	mmol/l	8.26	8.26	8.26	8.32	8.32	8.26	8.32	8.32	8.32	± 0.12
HT HCT	%	37.9	37.7	36.7	36.8	36.8	36.3	36.8	36.8	36.8	± 1.0
	l/l	0.379	0.377	0.367	0.368	0.368	0.363	0.368	0.368	0.368	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	242	249	240	242	245	237	245	244	242	± 10
VMP MPV	µm ³ ; fl	8.2	8.4	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300090308

ABX Minocal



LOT		CX 475		CAL		(Exp.) 2022-11-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE		
			YUMIZEN										
			H1500 H2500										
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	9.00									± 0.20	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58									± 0.06	
HB	HGB	g/dl	13.6									± 0.2	
		g/l	136									± 2	
		mmol/l	8.45									± 0.12	
HT	HCT	%	N/A									N/A	
		l/l	N/A									N/A	
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	244									± 10	
VMP	MPV	µm ³ ; fl	N/A									N/A	
VGM	MCV	µm ³ ; fl	83.0									± 4.0	

 1300060884
 Ref: TEMP-0387 Rev.48