

ABX Minocal


 (Exp.) 2023-01-05
 (YYYY - MM - DD)

LOT		CX 477				CAL		(Exp.) 2023-01-05 (YYYY - MM - DD)			
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.1	9.2	9.0	9.2	8.9	9.2	9.1	9.0	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.54	4.54	4.56	4.58	4.51	4.58	4.58	4.56	± 0.06
HB HGB	g/dl	13.9	13.8	13.8	14.0	14.0	13.9	14.0	14.0	14.0	± 0.2
	g/l	139	138	138	140	140	139	140	140	140	± 2
	mmol/l	8.63	8.57	8.57	8.69	8.69	8.63	8.69	8.69	8.69	± 0.12
HT HCT	%	39.4	39.0	38.4	38.3	38.2	37.9	38.2	38.7	38.3	± 1.0
	l/l	0.394	0.390	0.384	0.383	0.382	0.379	0.382	0.387	0.383	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	260	254	246	245	257	250	257	249	245	± 10
VMP MPV	µm ³ ; fl	8.9	8.9	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

 1300090308
 BACK / VERSO
 Ref. TEMP-0387 Rev.48