

ABX Minocal


 (Exp.) 2022-12-05
 (YYYY - MM - DD)

LOT		CX 476		CAL		(Exp.) 2022-12-05 (YYYY - MM - DD)					
Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.4	9.4	9.2	9.4	9.1	9.4	9.2	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.50	4.47	4.50	4.52	4.48	4.52	4.50	4.50	± 0.06
HB HGB	g/dl	13.2	13.2	13.1	13.2	13.2	13.2	13.2	13.2	13.2	± 0.2
	g/l	132	132	131	132	132	132	132	132	132	± 2
	mmol/l	8.20	8.20	8.14	8.20	8.20	8.20	8.20	8.20	8.20	± 0.12
HT HCT	%	37.2	37.4	36.2	36.2	36.2	36.1	36.2	36.2	36.2	± 1.0
	l/l	0.372	0.374	0.362	0.362	0.362	0.361	0.362	0.362	0.362	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	252	252	240	237	246	240	246	239	237	± 10
VMP MPV	µm ³ ; fl	9.1	9.2	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300090308