

ABX Minocal



LOT		CX 472		CAL		 (Exp.) 2022-08-05 (YYYY - MM - DD)								
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio							TOLERANCES TOLERANCE			
				MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA			PENTRA	PENTRA
				EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120			60 60 C+ ES60	XL80
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.2	9.3	9.0	9.3	9.2	9.3	9.0	9.0	± 0.2			
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.44	4.50	4.44	4.50	4.54	4.48	4.54	4.50	4.50	± 0.06			
HB HGB	g/dl	13.2	13.4	13.3	13.6	13.5	13.4	13.5	13.5	13.6	± 0.2			
	g/l	132	134	133	136	135	134	135	135	136	± 2			
	mmol/l	8.20	8.32	8.26	8.45	8.38	8.32	8.38	8.38	8.45	± 0.12			
HT HCT	%	37.7	38.3	36.6	37.4	37.2	36.7	37.2	36.9	37.4	± 1.0			
	l/l	0.377	0.383	0.366	0.374	0.372	0.367	0.372	0.369	0.374	± 0.010			
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	251	246	236	238	249	243	249	241	238	± 10			
VMP MPV	µm ³ ; fl	9.0	9.0	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5			

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300090308