

ABX Minocal



LOT

CX 472

Rev 1

CAL



(Exp.) 2022-08-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE	
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3					
		H550	H500 OT H500 CT			H500 OT	H500 CT H550				
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.90	8.90			9.00	9.00				± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57			4.48	4.48				± 0.06
HB HGB	g/dl	13.2	13.2			13.3	13.3				± 0.2
	g/l	132	132			133	133				± 2
	mmol/l	8.20	8.20			8.26	8.26				± 0.12
HT HCT	%	38.6	38.6			35.0	35.0				± 1.0
	l/l	0.386	0.386			0.350	0.350				± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	256	256			255	255				± 10
VMP MPV	μm ³ ; fl	10.9	10.9			10.6	10.6				± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208