

ABX Minocal



LOT

CX 474

CAL



(Exp.) 2022-10-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85	8.85			8.74	9.07			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.54			4.45	4.58			± 0.06
HB HGB	g/dl	12.9	12.9			12.8	13.0			± 0.2
	g/l	129	129			128	130			± 2
	mmol/l	8.01	8.01			7.95	8.07			± 0.12
HT HCT	%	37.7	37.7			36.2	36.6			± 1.0
	l/l	0.377	0.377			0.362	0.366			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	253	253			254	279			± 10
VMP MPV	μm ³ ; fl	10.4	10.4			10.7	10.8			± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208