

# ABX Minocal



**LOT**

**CX 473**

Rev 1

**CAL**



**(Exp.) 2022-09-05**  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE	
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3					
		H550	H500 OT H500 CT			H500 OT	H500 CT H550				
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.00</b>	<b>9.00</b>			<b>9.04</b>	<b>9.04</b>				± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.51</b>	<b>4.51</b>			<b>4.42</b>	<b>4.42</b>				± 0.06
HB HGB	g/dl	<b>12.8</b>	<b>12.8</b>			<b>12.8</b>	<b>12.8</b>				± 0.2
	g/l	<b>128</b>	<b>128</b>			<b>128</b>	<b>128</b>				± 2
	mmol/l	<b>7.95</b>	<b>7.95</b>			<b>7.97</b>	<b>7.97</b>				± 0.12
HT HCT	%	<b>37.4</b>	<b>37.4</b>			<b>35.8</b>	<b>35.8</b>				± 1.0
	l/l	<b>0.374</b>	<b>0.374</b>			<b>0.358</b>	<b>0.358</b>				± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>260</b>	<b>260</b>			<b>269</b>	<b>269</b>				± 10
VMP MPV	μm <sup>3</sup> ; fl	<b>10.3</b>	<b>10.3</b>			<b>10.4</b>	<b>10.4</b>				± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208