

ABX Minocal



LOT

CX 475

CAL



(Exp.) 2022-11-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.15	9.15			9.05	9.26			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.57			4.50	4.60			± 0.06
HB HGB	g/dl	13.1	13.1			13.1	13.2			± 0.2
	g/l	131	131			131	132			± 2
	mmol/l	8.14	8.14			8.14	8.20			± 0.12
HT HCT	%	38.2	38.2			36.5	36.5			± 1.0
	l/l	0.382	0.382			0.365	0.365			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	245	245			238	236			± 10
VMP MPV	µm ³ ; fl	10.4	10.4			10.3	10.3			± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208