

ABX Minocal



LOT

CX 476

CAL



(Exp.) 2022-12-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.93	8.93			8.56	9.26			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.54			4.45	4.57			± 0.06
HB HGB	g/dl	12.9	12.9			12.6	13.0			± 0.2
	g/l	129	129			126	130			± 2
	mmol/l	8.01	8.01			7.82	8.07			± 0.12
HT HCT	%	37.7	37.7			35.0	35.8			± 1.0
	l/l	0.377	0.377			0.350	0.358			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	249	249			259	265			± 10
VMP MPV	µm ³ ; fl	11.5	11.5			10.6	10.8			± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208