

ABX Minocal



LOT

CX 477

CAL



(Exp.) 2023-01-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.78	8.78			8.45	9.16			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.62			4.54	4.66			± 0.06
HB HGB	g/dl	13.6	13.6			13.4	13.7			± 0.2
	g/l	136	136			134	137			± 2
	mmol/l	8.45	8.45			8.32	8.51			± 0.12
HT HCT	%	39.3	39.3			37.3	38.3			± 1.0
	l/l	0.393	0.393			0.373	0.383			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	270	270			255	275			± 10
VMP MPV	μm ³ ; fl	11.1	11.1			10.6	11.0			± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208