



LOT		MX 437 L		Rev 1		CONTROL		L		 (Exp.)		2022-12-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.1	2.1			2.1						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.26	2.25			2.25						± 0.17	
HB HGB	g/dl	5.9	5.9			5.9						± 0.5	
	g/l	59	59			59						± 5	
	mmol/l	3.66	3.66			3.66						± 0.31	
HT HCT	%	15.8	15.8			15.8						± 2.0	
	l/l	0.158	0.158			0.158						± 0.020	
VGM MCV	µm ³ ; fl	70	70			70						± 4	
TGMH MCH	pg	26.1	26.2			26.2						± 2.0	
	fmol	1.62	1.63			1.63						± 0.12	
CCMH MCHC	g/dl	37.3	37.5			37.5						± 3.0	
	g/l	373	375			375						± 30	
	mmol/l	23.16	23.26			23.26						± 1.86	
IDR RDW	%	15.5	14.5			14.5						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	72	73			73						± 22	
VPM MPV	µm ³ ; fl	8.7	8.7			8.0						± 2.0	
LYMPHO	#	1.32	1.32			1.32						± 0.40	
	%	63.0	63.0			63.0						± 8.0	
MONO	#	0.27	0.27			0.27						± 0.20	
	%	13.0	13.0			13.0						± 6.0	
GRANULO	#	0.50	0.50			0.50						± 0.40	
	%	24.0	24.0			24.0						± 7.0	

LOT		MX 437 N		Rev 1		CONTROL N		(Exp.)		2022-12-05 (YYYY-MM-DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.7	7.6			7.7					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.60	4.60			4.59					± 0.26
HB HGB	g/dl	13.6	13.6			13.6					± 0.7
	g/l	136	136			136					± 7
	mmol/l	8.45	8.45			8.45					± 0.43
HT HCT	%	37.3	37.3			37.2					± 2.7
	l/l	0.373	0.373			0.372					± 0.027
VGM MCV	µm ³ ; fl	81	81			81					± 4
TGMH MCH	pg	29.6	29.6			29.6					± 2.0
	fmol	1.84	1.84			1.84					± 0.12
CCMH MCHC	g/dl	36.5	36.5			36.6					± 3.0
	g/l	365	365			366					± 30
	mmol/l	22.67	22.67			22.72					± 1.86
IDR RDW	%	15.0	14.0			14.0					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	269	272			261					± 40
VPM MPV	µm ³ ; fl	8.0	7.9			7.2					± 2.0
LYMPHO	#	2.58	2.70			2.62					± 0.40
	%	33.5	35.5			34.0					± 5.0
MONO	#	0.62	0.61			0.58					± 0.30
	%	8.0	8.0			7.5					± 4.0
GRANULO	#	4.50	4.29			4.50					± 0.60
	%	58.5	56.5			58.5					± 6.0

Ref: TEMP-0830 Rev.42 FRONT / RECTO 1300104144

LOT		MX 437 H		CONTROL		H		 (Exp.)		2022-12-05 (YYYY - MM - DD)				
		Rev 1		ABX Lyse										
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
				ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60						
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l		20.7	20.2			20.7						± 1.7
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.74	5.70			5.79						± 0.32
HB	HGB			18.3	18.2			18.3						± 0.9
		g/dl						18.3						± 9
		g/l					183						± 0.56	
		mmol/l		11.36	11.30			11.36						± 0.9
		%		51.1	51.3			52.1					± 3.4	
		l/l		0.511	0.513			0.521					± 0.034	
HT	HCT													
VGM	MCV	µm ³ ; fl		89	90			90						± 4
TGMH	MCH	pg		31.9	31.9			31.6						± 2.0
		fmol		1.98	1.98			1.96						± 0.12
CCMH	MCHC	g/dl		35.8	35.5			35.1						± 3.0
		g/l		358	355			351						± 30
		mmol/l		22.25	22.03			21.81						± 1.86
IDR	RDW	%		14.0	14.0			13.0						± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l		560	546			550						± 65
VPM	MPV	µm ³ ; fl		7.9	7.7			7.3						± 2.0
LYMPHO		#		3.31	3.23			3.31						± 1.00
		%		16.0	16.0			16.0						± 5.0
MONO		#		1.35	1.41			1.14						± 0.80
		%		6.5	7.0			5.5						± 4.0
GRANULO		#		16.04	15.55			16.25						± 1.40
		%		77.5	77.0			78.5						± 7.0