

LOT		MX 439 L		CONTROL		L		 (Exp.)		2023-04-05 (YYYY - MM - DD)													
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP 200								60 8P - 60 18P ADVIA 60									
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		2.1		1.9						2.0								± 0.4			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		2.33		2.35						2.31								± 0.17			
HB	HGB	g/dl		6.2		6.1						6.2								± 0.5			
		g/l		62		61						62								± 5			
		mmol/l		3.85		3.79						3.85								± 0.31			
HT	HCT	%		16.8		16.9						16.6								± 2.0			
		l/l		0.168		0.169						0.166								± 0.020			
VGM	MCV	µm ³ ; fl		72		72						72								± 4			
TGMH	MCH	pg		26.6		26.0						26.8								± 2.0			
		fmol		1.65		1.61						1.67								± 0.12			
CCMH	MCHC	g/dl		37.0		36.1						37.3								± 3.0			
		g/l		370		361						373								± 30			
		mmol/l		22.95		22.39						23.15								± 1.86			
IDR	RDW	%		15.5		14.0						14.0								± 3.0			
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		78		79						76								± 22			
VPM	MPV	µm ³ ; fl		9.2		8.5						8.4								± 2.0			
LYMPHO		#		1.34		1.24						1.28								± 0.40			
		%		64.0		65.0						64.0								± 8.0			
MONO		#		0.19		0.17						0.18								± 0.20			
		%		9.0		9.0						9.0								± 6.0			
GRANULO		#		0.57		0.49						0.54								± 0.40			
		%		27.0		26.0						27.0								± 7.0			

LOT		MX 439 N		CONTROL		N		(Exp.)		2023-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.6	7.3			7.6					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.53	4.52			4.49					± 0.26
HB HGB	g/dl	13.4	13.5			13.4					± 0.7
	g/l	134	135			134					± 7
	mmol/l	8.32	8.38			8.32					± 0.43
HT HCT	%	36.7	37.1			36.8					± 2.7
	l/l	0.367	0.371			0.368					± 0.027
VGM MCV	µm ³ ; fl	81	82			82					± 4
TGMH MCH	pg	29.6	29.9			29.8					± 2.0
	fmol	1.84	1.85			1.85					± 0.12
CCMH MCHC	g/dl	36.5	36.4			36.4					± 3.0
	g/l	365	364			364					± 30
	mmol/l	22.68	22.62			22.60					± 1.86
IDR RDW	%	15.0	13.5			13.0					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	278	272			268					± 40
VPM MPV	µm ³ ; fl	8.6	8.2			7.9					± 2.0
LYMPHO	#	2.51	2.45			2.51					± 0.40
	%	33.0	33.5			33.0					± 5.0
MONO	#	0.57	0.62			0.53					± 0.30
	%	7.5	8.5			7.0					± 4.0
GRANULO	#	4.52	4.23			4.56					± 0.60
	%	59.5	58.0			60.0					± 6.0

Ref: TEMP-0830 Rev.42 FRONT / RECTO 1300104144



LOT		MX 439 H		CONTROL		H		 (Exp.)		2023-04-05 (YYYY - MM - DD)								
		Rev 1		ABX Lyse														
PARAMETRES PARAMETERS	UNITES UNITS											TOLERANCES TOLERANCE						
		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS			PENTRA		PENTRA		PENTRA	
		ES60 Care ST		CRP 200				60 8P - 60 18P ADVIA 60										
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.5		19.6				20.3								± 1.7		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.54		5.50				5.50								± 0.32		
HB HGB	g/dl	17.7		17.6				17.7								± 0.9		
	g/l	177		176				177								± 9		
	mmol/l	10.99		10.93				10.99								± 0.56		
HT HCT	%	49.3		49.0				49.5								± 3.4		
	l/l	0.493		0.490				0.495								± 0.034		
VGM MCV	µm ³ ; fl	89		89				90								± 4		
TGMH MCH	pg	31.9		32.0				32.2								± 2.0		
	fmol	1.98		1.99				2.00								± 0.12		
CCMH MCHC	g/dl	35.9		36.0				35.8								± 3.0		
	g/l	359		360				358								± 30		
	mmol/l	22.29		22.33				22.21								± 1.86		
IDR RDW	%	14.0		13.0				12.5								± 3.0		
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	530		525				516								± 65		
VPM MPV	µm ³ ; fl	8.5		8.2				7.8								± 2.0		
LYMPHO	#	3.49		3.23				3.25								± 1.00		
	%	17.0		16.5				16.0								± 5.0		
MONO	#	1.13		1.18				1.12								± 0.80		
	%	5.5		6.0				5.5								± 4.0		
GRANULO	#	15.89		15.19				15.94								± 1.40		
	%	77.5		77.5				78.5								± 7.0		