

ABX Minocal



LOT

CX 478

CAL



(Exp.) 2023-02-05

Rev **2**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.0			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.60			4.55					± 0.06
HB HGB	g/dl	13.5	13.7			13.5					± 0.2
	g/l	135	137			135					± 2
	mmol/l	8.38	8.51			8.38					± 0.12
HT HCT	%	38.5	38.9			38.4					± 1.0
	l/l	0.385	0.389			0.384					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	242	242			245					± 10
VMP MPV	µm ³ ; fl	8.7	8.3			7.9					± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300090308