

ABX Minocal



LOT		CX 478		CAL		(Exp.) 2023-02-05 (YYYY - MM - DD)					
		Rev 2									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.0			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.60			4.53					± 0.06
HB HGB	g/dl	13.6	13.4			13.7					± 0.2
	g/l	136	134			137					± 2
	mmol/l	8.45	8.34			8.51					± 0.12
HT HCT	%	38.1	38.9			39.2					± 1.0
	l/l	0.381	0.389			0.392					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	227	242			241					± 10
VMP MPV	µm ³ ; fl	8.7	8.3			9.1					± 0.5

ABX Minocal



LOT		CX 478		CAL		(Exp.) 2023-02-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.8								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.53								± 0.06
HB HGB	g/dl	13.2	13.4								± 0.2
	g/l	132	134								± 2
	mmol/l	8.20	8.32								± 0.12
HT HCT	%	37.7	38.5								± 1.0
	l/l	0.377	0.385								± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	219	245								± 10
VMP MPV	µm ³ ; fl	N/A	9.2								± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208