

ABX Minocal



(Exp.) 2023-02-05

(YYYY - MM - DD)

LOT		CX 478		CAL								
		Rev 2				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.0			9.1					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.60			4.55					± 0.06	
HB HGB	g/dl	13.5	13.7			13.5					± 0.2	
	g/l	135	137			135					± 2	
	mmol/l	8.38	8.51			8.38					± 0.12	
HT HCT	%	38.5	38.9			38.4					± 1.0	
	l/l	0.385	0.389			0.384					± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	242	242			245					± 10	
VMP MPV	μm ³ ; fl	8.7	8.3			7.9					± 0.5	

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.48