

ABX Minocal


 (Exp.) **2023-02-05**
 (YYYY - MM - DD)

LOT		CX 478				CAL		(Exp.) 2023-02-05 (YYYY - MM - DD)			
		Rev 2									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60	DX120 - DF120	60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2	9.3	9.1	9.3	9.1	9.3	9.1	9.1	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.53	4.47	4.52	4.54	4.51	4.54	4.51	4.52	± 0.06
HB HGB	g/dl	13.4	13.4	13.4	13.5	13.5	13.5	13.5	13.5	13.5	± 0.2
	g/l	134	134	134	135	135	135	135	135	135	± 2
	mmol/l	8.32	8.32	8.32	8.38	8.38	8.38	8.38	8.38	8.38	± 0.12
HT HCT	%	37.8	37.8	36.2	36.8	36.8	36.5	36.8	37.2	36.8	± 1.0
	l/l	0.378	0.378	0.362	0.368	0.368	0.365	0.368	0.372	0.368	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	250	244	235	239	246	238	246	240	239	± 10
VMP MPV	µm ³ ; fl	8.7	8.7	N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300090308