

# ABX Minocal


 (Exp.) 2023-02-05  
 (YYYY - MM - DD)

| LOT                      |                                                       | CX 478      |         | CAL    |        | <div><div></div><div>(Exp.) 2023-02-05</div><div>(YYYY - MM - DD )</div></div> |        |               |                     |        |                         |
|--------------------------|-------------------------------------------------------|-------------|---------|--------|--------|--------------------------------------------------------------------------------|--------|---------------|---------------------|--------|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio |         |        |        |                                                                                |        |               |                     |        | TOLERANCES<br>TOLERANCE |
|                          |                                                       | MICROS      | MICROS  | PENTRA | PENTRA | PENTRA                                                                         | PENTRA | PENTRA        | PENTRA              | PENTRA |                         |
|                          |                                                       | EMI CRP     | LC-767G | MS CRP | XLR    | DX NEXUS<br>DF NEXUS                                                           | MS60   | DX120 - DF120 | 60<br>60 C+<br>ES60 | XL80   |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.2         | 9.2     | 9.3    | 9.1    | 9.3                                                                            | 9.1    | 9.3           | 9.1                 | 9.1    | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.53        | 4.53    | 4.47   | 4.52   | 4.54                                                                           | 4.51   | 4.54          | 4.51                | 4.52   | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.4        | 13.4    | 13.4   | 13.5   | 13.5                                                                           | 13.5   | 13.5          | 13.5                | 13.5   | ± 0.2                   |
|                          | g/l                                                   | 134         | 134     | 134    | 135    | 135                                                                            | 135    | 135           | 135                 | 135    | ± 2                     |
|                          | mmol/l                                                | 8.32        | 8.32    | 8.32   | 8.38   | 8.38                                                                           | 8.38   | 8.38          | 8.38                | 8.38   | ± 0.12                  |
| HT HCT                   | %                                                     | 37.8        | 37.8    | 36.2   | 36.8   | 36.8                                                                           | 36.5   | 36.8          | 37.2                | 36.8   | ± 1.0                   |
|                          | l/l                                                   | 0.378       | 0.378   | 0.362  | 0.368  | 0.368                                                                          | 0.365  | 0.368         | 0.372               | 0.368  | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 250         | 244     | 235    | 239    | 246                                                                            | 238    | 246           | 240                 | 239    | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 8.7         | 8.7     | N/A    | N/A    | N/A                                                                            | N/A    | N/A           | N/A                 | N/A    | ± 0.5                   |

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300090308