

ABX Minocal



(Exp.) 2023-03-05

(YYYY - MM - DD)

LOT		CX 479		CAL		(Exp.) 2023-03-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.55			4.52					± 0.06
HB HGB	g/dl	13.4	13.6			13.5					± 0.2
	g/l	134	136			135					± 2
	mmol/l	8.32	8.45			8.38					± 0.12
HT HCT	%	38.2	38.2			38.4					± 1.0
	l/l	0.382	0.382			0.384					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	256	250			250					± 10
VMP MPV	μm ³ ; fl	8.7	8.9			8.2					± 0.5

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.48