

ABX Minocal

LOT		CX 479		CAL		(Exp.) 2023-03-05 (YYYY - MM - DD)						
		Rev 2										
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA		
		ES60	CRP 200			60 8P 60 18P						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9			8.9					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.55			4.55					± 0.06	
HB HGB	g/dl	13.5	13.3			13.4					± 0.2	
	g/l	135	133			134					± 2	
	mmol/l	8.38	8.28			8.32					± 0.12	
HT HCT	%	37.9	38.2			38.7					± 1.0	
	l/l	0.379	0.382			0.387					± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234	250			252					± 10	
VMP MPV	µm ³ ; fl	9.6	8.9			9.3					± 0.5	

ABX Minocal



LOT		CX 479		CAL		(Exp.) 2023-03-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	LyseVet									TOLERANCES TOLERANCE
		MICROS	MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7	8.7								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.55								± 0.06
HB HGB	g/dl	13.0	13.2								± 0.2
	g/l	130	132								± 2
	mmol/l	8.07	8.20								± 0.12
HT HCT	%	37.4	38.2								± 1.0
	l/l	0.374	0.382								± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	229	245								± 10
VMP MPV	µm ³ ; fl	N/A	9.4								± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208