

ABX Minocal

LOT

CX 479

CAL



(Exp.) 2023-03-05

Rev **2**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.55			4.55					± 0.06
HB HGB	g/dl	13.5	13.3			13.4					± 0.2
	g/l	135	133			134					± 2
	mmol/l	8.38	8.28			8.32					± 0.12
HT HCT	%	37.9	38.2			38.7					± 1.0
	l/l	0.379	0.382			0.387					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234	250			252					± 10
VMP MPV	µm ³ ; fl	9.6	8.9			9.3					± 0.5

ABX Minocal

LOT

CX 479

CAL



(Exp.) 2023-03-05
(YYYY-MM-DD)

Rev **2**

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.55							± 0.06
HB HGB	g/dl	13.0	13.2							± 0.2
	g/l	130	132							± 2
	mmol/l	8.07	8.20							± 0.12
HT HCT	%	37.4	38.2							± 1.0
	l/l	0.374	0.382							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	229	245							± 10
VMP MPV	µm ³ ; fl	N/A	9.4							± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208