

ABX Minocal



LOT		CX 480		CAL		(Exp.) 2023-04-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7	8.8			8.7					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.59			4.55					± 0.06
HB HGB	g/dl	13.5	13.3			13.7					± 0.2
	g/l	135	133			137					± 2
	mmol/l	8.38	8.28			8.51					± 0.12
HT HCT	%	38.3	38.8			39.8					± 1.0
	l/l	0.383	0.388			0.398					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233	246			249					± 10
VMP MPV	µm ³ ; fl	9.1	8.6			9.2					± 0.5

ABX Minocal



LOT		CX 480		CAL		(Exp.) 2023-04-05 (YYYY - MM - DD)						
		Rev 1	LyseVet									
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS							TOLERANCES TOLERANCE
				Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC		10 ³ /mm ³ ;10 ⁹ /l	8.6	8.5								± 0.2
GR RBC		10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.52								± 0.06
HB HGB		g/dl	13.3	13.4								± 0.2
		g/l	133	134								± 2
		mmol/l	8.26	8.32								± 0.12
HT HCT		%	38.2	38.4								± 1.0
		l/l	0.382	0.384								± 0.010
PLA PLT		10 ³ /mm ³ ;10 ⁹ /l	220	253								± 10
VMP MPV		µm ³ ; fl	N/A	8.9								± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208