

ABX Minocal



(Exp.) 2023-04-05

(YYYY - MM - DD)

LOT		CX 480		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7	8.8			8.7					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.59			4.58					± 0.06	
HB HGB	g/dl	13.5	13.6			13.6					± 0.2	
	g/l	135	136			136					± 2	
	mmol/l	8.38	8.45			8.45					± 0.12	
HT HCT	%	38.5	38.8			38.7					± 1.0	
	l/l	0.385	0.388			0.387					± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	248	246			245					± 10	
VMP MPV	µm ³ ; fl	8.7	8.6			8.2					± 0.5	

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.48