

ABX Minocal

LOT
CX 481
CAL

(Exp.) 2023-05-05

(YYYY - MM - DD)

 Rev **2**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.8			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39	4.38			4.39					± 0.06
HB HGB	g/dl	12.9	12.9			13.0					± 0.2
	g/l	129	129			130					± 2
	mmol/l	8.01	8.01			8.07					± 0.12
HT HCT	%	36.4	36.4			36.7					± 1.0
	l/l	0.364	0.364			0.367					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	244	244			243					± 10
VMP MPV	µm ³ ; fl	8.7	8.6			8.0					± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300090308