

# ABX Minotrol CRP



**LOT**

**MC 427**

**CONTROL**



**(Exp.) 2023-07-05**  
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS		UNITES UNITS		Minilysebio												TOLERANCES TOLERANCE	
				CONTROL			1	CONTROL			2	CONTROL			3		
				MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE	MICROS			TOLERANCES TOLERANCE		
				CRP 200	60 8P 60 18P			CRP 200	60 8P 60 18P			CRP 200	60 8P 60 18P				
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	2.0	2.0		± 0.4	7.9	8.0		± 1.0	19.6	20.2		± 2.2				
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	2.34	2.33		± 0.15	4.47	4.46		± 0.18	5.67	5.60		± 0.25				
HB HGB	g/dl	5.8	5.8		± 0.4	13.1	13.3		± 0.5	18.3	18.8		± 0.7				
	g/l	58	58		± 4	131	133		± 5	183	188		± 7				
	mmol/l	3.60	3.60		± 0.25	8.15	8.26		± 0.31	11.36	11.67		± 0.44				
HT HCT	%	16.6	16.5		± 2.0	36.7	37.5		± 2.5	51.6	52.6		± 3.0				
	l/l	0.166	0.165		± 0.020	0.367	0.375		± 0.025	0.516	0.526		± 0.030				
VGM MCV	µm <sup>3</sup> ; fl	71	71		± 4	82	84		± 4	91	94		± 4				
TGMH MCH	pg	24.8	24.9		± 2.4	29.4	29.8		± 2.4	32.3	33.6		± 2.8				
	fmol	1.54	1.55		± 0.15	1.82	1.85		± 0.15	2.00	2.08		± 0.18				
CCMH MCHC	g/dl	34.9	35.1		± 3.0	35.8	35.5		± 3.0	35.5	35.7		± 3.0				
	g/l	349	351		± 30	358	355		± 30	355	357		± 30				
	mmol/l	21.66	21.77		± 1.86	22.25	22.05		± 1.86	22.03	22.18		± 1.86				
IDR RDW	%	14.0	14.5		± 3.0	13.5	14.0		± 3.0	13.0	13.5		± 3.0				
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	76	72		± 20	256	250		± 40	509	515		± 60				
VPM MPV	µm <sup>3</sup> ; fl	8.7	8.9		± 2.0	8.2	8.5		± 2.0	7.9	8.3		± 2.0				
LYMPHO	#	1.4	1.4		± 0.2	2.8	2.7		± 0.5	3.6	3.5		± 1.1				
	%	70.7	68.0		± 8.0	35.2	33.5		± 6.0	18.2	17.5		± 5.0				
MONO	#	0.2	0.2		± 0.2	0.6	0.4		± 0.4	1.1	0.8		± 1.0				
	%	7.6	7.5		± 7.0	7.6	5.5		± 5.0	5.5	4.0		± 5.0				
GRANULO	#	0.4	0.5		± 0.2	4.5	4.9		± 0.5	15.0	15.9		± 1.5				
	%	22.1	24.5		± 7.0	57.3	61.0		± 6.0	76.7	78.5		± 7.0				
CRP	mg/dl	0.30	N/A		± 0.30	2.08	N/A		± 0.70	4.08	N/A		± 1.00				
	mg/l	3.0	N/A		± 3.0	20.8	N/A		± 7.0	40.8	N/A		± 10.0				

1300027267-A  
Ref: TEMP-1070 Rev.42