

ABX Minocal

LOT

CX 481

CAL



(Exp.) 2023-05-05

Rev **2**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.8			8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39	4.38			4.37					± 0.06
HB HGB	g/dl	13.0	12.6			12.9					± 0.2
	g/l	130	126			129					± 2
	mmol/l	8.07	7.85			8.01					± 0.12
HT HCT	%	36.4	36.4			36.9					± 1.0
	l/l	0.364	0.364			0.369					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	225	244			246					± 10
VMP MPV	µm ³ ; fl	9.3	8.6			9.2					± 0.5

ABX Minocal

LOT

CX 481

CAL



(Exp.) 2023-05-05
(YYYY-MM-DD)

Rev **2**

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.8	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.40							± 0.06
HB HGB	g/dl	12.7	12.7							± 0.2
	g/l	127	127							± 2
	mmol/l	7.89	7.89							± 0.12
HT HCT	%	36.3	36.5							± 1.0
	l/l	0.363	0.365							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	222	237							± 10
VMP MPV	µm ³ ; fl	N/A	8.9							± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208