

# ABX Minocal


 (Exp.) 2023-05-05  
 (YYYY - MM - DD)

| LOT                      |                                                       | CX 481      |         | CAL    |        | <div><div><div></div></div><div>(Exp.) 2023-05-05</div><div>(YYYY - MM - DD )</div></div> |        |               |                     |        |                         |
|--------------------------|-------------------------------------------------------|-------------|---------|--------|--------|-------------------------------------------------------------------------------------------|--------|---------------|---------------------|--------|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio |         |        |        |                                                                                           |        |               |                     |        | TOLERANCES<br>TOLERANCE |
|                          |                                                       | MICROS      | MICROS  | PENTRA | PENTRA | PENTRA                                                                                    | PENTRA | PENTRA        | PENTRA              | PENTRA |                         |
|                          |                                                       | EMI CRP     | LC-767G | MS CRP | XLR    | DX NEXUS<br>DF NEXUS                                                                      | MS60   | DX120 - DF120 | 60<br>60 C+<br>ES60 | XL80   |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 8.9         | 9.0     | 9.1    | 9.0    | 9.1                                                                                       | 8.9    | 9.1           | 9.1                 | 9.0    | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.33        | 4.37    | 4.39   | 4.38   | 4.40                                                                                      | 4.35   | 4.40          | 4.41                | 4.38   | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 12.7        | 12.7    | 12.8   | 12.9   | 12.8                                                                                      | 12.8   | 12.8          | 12.9                | 12.9   | ± 0.2                   |
|                          | g/l                                                   | 127         | 127     | 128    | 129    | 128                                                                                       | 128    | 128           | 129                 | 129    | ± 2                     |
|                          | mmol/l                                                | 7.89        | 7.89    | 7.95   | 8.01   | 7.95                                                                                      | 7.95   | 7.95          | 8.01                | 8.01   | ± 0.12                  |
| HT HCT                   | %                                                     | 35.5        | 36.1    | 35.6   | 35.5   | 35.0                                                                                      | 34.8   | 35.0          | 35.5                | 35.5   | ± 1.0                   |
|                          | l/l                                                   | 0.355       | 0.361   | 0.356  | 0.355  | 0.350                                                                                     | 0.348  | 0.350         | 0.355               | 0.355  | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 244         | 246     | 239    | 236    | 250                                                                                       | 240    | 250           | 245                 | 236    | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 8.7         | 8.8     | N/A    | N/A    | N/A                                                                                       | N/A    | N/A           | N/A                 | N/A    | ± 0.5                   |

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300090308