

ABX Minocal



LOT

CX 481
Rev 2

CAL

(Exp.) **2023-05-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80			8.80	8.80			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42			4.45	4.45			± 0.06
HB HGB	g/dl	12.5	12.5			12.5	12.5			± 0.2
	g/l	125	125			125	125			± 2
	mmol/l	7.76	7.76			7.76	7.76			± 0.12
HT HCT	%	36.5	36.5			35.2	34.4			± 1.0
	l/l	0.365	0.365			0.352	0.344			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	248	248			275	278			± 10
VMP MPV	μm ³ ; fl	10.7	10.7			9.2	9.2			± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208