



Yumizen G PT Liq 4 for Yumizen G Line



Reference No.:		1300036373		
Lot No.:		920315 (rev. 1)		
Expiry date:		2023-09-30		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		%	sec
		100	11,9	
		50	16,3	
		33	20,7	
		25	24,9	
		1,19	MNPT 11,9	
Important notice				
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p>				