

Yumizen G PT Liq 4 for Yumizen G Line



| | | | | |
|--|-----------------------|------------------------|---------------------|-------------|
| Reference No.: | | 1300036373 | | |
| Lot No.: | | 920611 (rev. 1) | | |
| Expiry date: | | 2023-12-31 | | |
| Coagulometers | | ISI | MASTER CURVE | |
| Optical | Yumizen G Line | | % | sec |
| | | 100 | 13,3 | |
| | | 50 | 19,4 | |
| | | 33 | 25,7 | |
| | | 25 | 31,6 | |
| | | 1,15 | MNPT | 13,3 |
| Important notice | | | | |
| <p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p> | | | | |