

LOT		MX 441 L		Rev 1		CONTROL		L		 (Exp.)		2023-08-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	2.0	2.0			2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.31	2.32			2.29						± 0.17	
HB HGB	g/dl	6.1	6.1			6.1						± 0.5	
	g/l	61	61			61						± 5	
	mmol/l	3.79	3.79			3.79						± 0.31	
HT HCT	%	16.4	16.2			16.3						± 2.0	
	l/l	0.164	0.162			0.163						± 0.020	
VGM MCV	µm ³ ; fl	71	70			71						± 4	
TGMH MCH	pg	26.4	26.3			26.6						± 2.0	
	fmol	1.64	1.63			1.65						± 0.12	
CCMH MCHC	g/dl	37.2	37.6			37.5						± 3.0	
	g/l	372	376			375						± 30	
	mmol/l	23.10	23.33			23.30						± 1.86	
IDR RDW	%	14.5	13.5			13.5						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	77	85			80						± 22	
VPM MPV	µm ³ ; fl	8.5	8.0			7.8						± 2.0	
LYMPHO	#	1.28	1.33			1.30						± 0.40	
	%	64.0	66.5			65.0						± 8.0	
MONO	#	0.19	0.19			0.19						± 0.20	
	%	9.5	9.5			9.5						± 6.0	
GRANULO	#	0.53	0.48			0.51						± 0.40	
	%	26.5	24.0			25.5						± 7.0	

LOT		MX 441 N		CONTROL		N		 (Exp.)		2023-08-05 (YYYY - MM - DD)											
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP 200				60 8P - 60 18P ADVIA 60											
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.5	7.4				7.5											± 0.8			
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.49	4.52				4.51											± 0.26			
HB HGB	g/dl	13.5	13.4				13.5											± 0.7			
	g/l	135	134				135											± 7			
	mmol/l	8.38	8.32				8.38											± 0.43			
HT HCT	%	36.8	36.6				37.0											± 2.7			
	l/l	0.368	0.366				0.370											± 0.027			
VGM MCV	µm ³ ; fl	82	81				82											± 4			
TGMH MCH	pg	30.1	29.6				29.9											± 2.0			
	fmol	1.87	1.84				1.86											± 0.12			
CCMH MCHC	g/dl	36.7	36.6				36.5											± 3.0			
	g/l	367	366				365											± 30			
	mmol/l	22.77	22.73				22.67											± 1.86			
IDR RDW	%	14.5	13.5				13.0											± 3.0			
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	265	254				260											± 40			
VPM MPV	µm ³ ; fl	8.2	8.0				7.4											± 2.0			
LYMPHO	#	2.51	2.59				2.59											± 0.40			
	%	33.5	35.0				34.5											± 5.0			
MONO	#	0.56	0.59				0.56											± 0.30			
	%	7.5	8.0				7.5											± 4.0			
GRANULO	#	4.43	4.22				4.35											± 0.60			
	%	59.0	57.0				58.0											± 6.0			

Rel: TEMP-0830 Rev.42 FRONT / RECTO 1300104144



LOT		MX 441 H		CONTROL		H		 (Exp.)		2023-08-05 (YYYY - MM - DD)											
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse							TOLERANCES TOLERANCE										
				MICROS		MICROS		MICROS		MICROS			ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		
				ES60 Care ST		CRP 200							60 8P - 60 18P ADVIA 60								
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.3	19.1					19.6									± 1.7			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.54	5.60					5.59									± 0.32			
HB	HGB	g/dl	17.7	17.8					17.9									± 0.9			
		g/l	177	178					179								± 9				
		mmol/l	10.99	11.05					11.12								± 0.56				
HT	HCT	%	49.9	49.8					50.3								± 3.4				
		l/l	0.499	0.498					0.503								± 0.034				
VGM	MCV	µm ³ ; fl	90	89					90								± 4				
TGMH	MCH	pg	31.9	31.8					32.0								± 2.0				
		fmol	1.98	1.97					1.99								± 0.12				
CCMH	MCHC	g/dl	35.5	35.7					35.6								± 3.0				
		g/l	355	357					356								± 30				
		mmol/l	22.05	22.18					22.09								± 1.86				
IDR	RDW	%	14.0	13.5					13.0								± 3.0				
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	526	507					520								± 65				
VPM	MPV	µm ³ ; fl	8.0	7.7					7.2								± 2.0				
LYMPHO	#	3.09	3.25					3.14									± 1.00				
	%	16.0	17.0					16.0									± 5.0				
MONO	#	0.97	1.15					0.98									± 0.80				
	%	5.0	6.0					5.0									± 4.0				
GRANULO	#	15.25	14.71					15.48									± 1.40				
	%	79.0	77.0					79.0									± 7.0				