

ABX Minocal



(Exp.) 2023-06-05

(YYYY - MM - DD)

LOT		CX 482		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST	CRP 200			60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.0			9.2						± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39	4.42			4.39						± 0.06
HB HGB	g/dl	13.1	13.2			13.2						± 0.2
	g/l	131	132			132						± 2
	mmol/l	8.14	8.20			8.20						± 0.12
HT HCT	%	37.1	37.3			37.5						± 1.0
	l/l	0.371	0.373			0.375						± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	247	248			243						± 10
VMP MPV	µm ³ ; fl	8.5	8.2			7.6						± 0.5

 1300090308
 FRONT / RECTO
 Ref. TEMP-0387 Rev.48