

ABX Minocal

LOT		CX 482		CAL		(Exp.) 2023-06-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA		
		ES60	CRP 200			60 8P 60 18P						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0			9.2					± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39	4.42			4.36					± 0.06	
HB HGB	g/dl	13.2	12.9			13.2					± 0.2	
	g/l	132	129			132					± 2	
	mmol/l	8.20	8.03			8.20					± 0.12	
HT HCT	%	37.1	37.3			37.5					± 1.0	
	l/l	0.371	0.373			0.375					± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	229	248			251					± 10	
VMP MPV	µm ³ ; fl	8.5	8.2			8.8					± 0.5	

ABX Minocal



LOT		CX 482		CAL		(Exp.) 2023-06-05 (YYYY - MM - DD)						
		Rev 1	LyseVet									
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS							TOLERANCES TOLERANCE
				Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC		10 ³ /mm ³ ;10 ⁹ /l	9.2	8.9								± 0.2
GR RBC		10 ⁶ /mm ³ ;10 ¹² /l	4.38	4.40								± 0.06
HB HGB		g/dl	13.0	13.0								± 0.2
		g/l	130	130								± 2
		mmol/l	8.07	8.07								± 0.12
HT HCT		%	36.8	37.2								± 1.0
		l/l	0.368	0.372								± 0.010
PLA PLT		10 ³ /mm ³ ;10 ⁹ /l	227	248								± 10
VMP MPV		µm ³ ; fl	N/A	8.7								± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208