

ABX Minocal

LOT

CX 482

CAL



(Exp.) 2023-06-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0			9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.39	4.42			4.36					± 0.06
HB HGB	g/dl	13.2	12.9			13.2					± 0.2
	g/l	132	129			132					± 2
	mmol/l	8.20	8.03			8.20					± 0.12
HT HCT	%	37.1	37.3			37.5					± 1.0
	l/l	0.371	0.373			0.375					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	229	248			251					± 10
VMP MPV	µm ³ ; fl	8.5	8.2			8.8					± 0.5

ABX Minocal

LOT

CX 482

CAL



(Exp.) 2023-06-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.2	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.38	4.40							± 0.06
HB HGB	g/dl	13.0	13.0							± 0.2
	g/l	130	130							± 2
	mmol/l	8.07	8.07							± 0.12
HT HCT	%	36.8	37.2							± 1.0
	l/l	0.368	0.372							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	227	248							± 10
VMP MPV	µm ³ ; fl	N/A	8.7							± 0.5

Ref: TEMP-0387 Rev.48 FRONT / RECTO 1300105208