

# ABX Minocal


 (Exp.) 2023-06-05  
 (YYYY - MM - DD)

| LOT                      |                                                       | CX 482      |         |        |        | CAL                  |        | <div><div></div><div>(Exp.) 2023-06-05</div><div>(YYYY - MM - DD )</div></div> |                     |        |                         |
|--------------------------|-------------------------------------------------------|-------------|---------|--------|--------|----------------------|--------|--------------------------------------------------------------------------------|---------------------|--------|-------------------------|
|                          |                                                       | Rev 1       |         |        |        |                      |        |                                                                                |                     |        |                         |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio |         |        |        |                      |        |                                                                                |                     |        | TOLERANCES<br>TOLERANCE |
|                          |                                                       | MICROS      | MICROS  | PENTRA | PENTRA | PENTRA               | PENTRA | PENTRA                                                                         | PENTRA              | PENTRA |                         |
|                          |                                                       | EMI CRP     | LC-767G | MS CRP | XLR    | DX NEXUS<br>DF NEXUS | MS60   | DX120 - DF120                                                                  | 60<br>60 C+<br>ES60 | XL80   |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.2         | 9.4     | 9.3    | 9.1    | 9.4                  | 9.2    | 9.4                                                                            | 9.2                 | 9.1    | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.39        | 4.38    | 4.40   | 4.35   | 4.41                 | 4.34   | 4.41                                                                           | 4.38                | 4.35   | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.1        | 13.0    | 13.1   | 13.2   | 13.1                 | 13.1   | 13.1                                                                           | 13.2                | 13.2   | ± 0.2                   |
|                          | g/l                                                   | 131         | 130     | 131    | 132    | 131                  | 131    | 131                                                                            | 132                 | 132    | ± 2                     |
|                          | mmol/l                                                | 8.14        | 8.07    | 8.14   | 8.20   | 8.14                 | 8.14   | 8.14                                                                           | 8.20                | 8.20   | ± 0.12                  |
| HT HCT                   | %                                                     | 36.9        | 36.8    | 35.9   | 36.1   | 35.9                 | 35.6   | 35.9                                                                           | 35.9                | 36.1   | ± 1.0                   |
|                          | l/l                                                   | 0.369       | 0.368   | 0.359  | 0.361  | 0.359                | 0.356  | 0.359                                                                          | 0.359               | 0.361  | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 253         | 247     | 238    | 243    | 253                  | 249    | 253                                                                            | 249                 | 243    | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 8.5         | 8.5     | N/A    | N/A    | N/A                  | N/A    | N/A                                                                            | N/A                 | N/A    | ± 0.5                   |

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300090308