

ABX Minocal



LOT

CX 482

CAL



(Exp.) 2023-06-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN V1.0 to V2.x				YUMIZEN since V3				
		H550	H500 OT H500 CT			H500 OT	H500 CT H550			
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.00	9.00			9.00	9.00			± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42			4.43	4.43			± 0.06
HB HGB	g/dl	12.8	12.8			12.8	12.8			± 0.2
	g/l	128	128			128	128			± 2
	mmol/l	7.95	7.95			7.95	7.95			± 0.12
HT HCT	%	37.1	37.1			36.1	35.1			± 1.0
	l/l	0.371	0.371			0.361	0.351			± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	255	255			275	280			± 10
VMP MPV	µm ³ ; fl	10.4	10.4			9.2	9.2			± 0.5

Ref: TEMP-0387 Rev.48 BACK / VERSO 1300105208