



**Yumizen G CTRL DDi I
 for YUMIZEN
 G1500/G1550**



D 4 9 2 0 9 3 1 4 6 4 1 4 0 9 2 4 D



D 5 0 5 2 1 0 2 0 7 7 D

Reference No.:		1300036414
Lot No.:		920931
Expiry date:		2024-09-30
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,52-1,02

**Yumizen G CTRL DDi I
 for YUMIZEN G800**



D 4 9 2 0 9 3 1 4 6 4 1 4 0 9 2 4 D



D 5 0 5 5 1 0 5 0 8 3 D

Reference No.:		1300036414
Lot No.:		920931
Expiry date:		2024-09-30
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,55-1,05

**Yumizen G CTRL DDi I
 for YUMIZEN
 G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		920931
Expiry date:		2024-09-30
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,46-0,96

***Yumizen G CTRL DDi II
 for YUMIZEN
 G1500/G1550***



D 4 9 2 0 9 3 1 5 6 4 1 4 0 9 2 4 D



D 5 1 3 8 2 4 9 0 9 5 D

Reference No.:		1300036414
Lot No.:		920931
Expiry date:		2024-09-30
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,38-2,49

***Yumizen G CTRL DDi II
 for YUMIZEN G800***



D 4 9 2 0 9 3 1 5 6 4 1 4 0 9 2 4 D



D 5 1 3 0 2 3 5 0 8 2 D

Reference No.:		1300036414
Lot No.:		920931
Expiry date:		2024-09-30
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,30-2,35

***Yumizen G CTRL DDi II
 for YUMIZEN
 G400 DDi/G200***

Reference No.:		1300036414
Lot No.:		920931
Expiry date:		2024-09-30
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,45-2,60