



Yumizen G FIB 2 for Yumizen G1500/G1550 and Yumizen G800

Reference No.:		1300036383	
Lot No.:		920704 (rev. 1)	
Expiry date:		2024-07-31	
Coagulometers		MASTER CURVE	
		g/L	sec
Optical	Yumizen G1500/G1550/G800	4	8,3
		3	10,0
		2	12,9
		1	19,9
		Dilution	10x



Yumizen G FIB 2 for Yumizen G400 DDi/G200

Reference No.:		1300036383	
Lot No.:		920704 (rev. 1)	
Expiry date:		2024-07-31	
Coagulometers		MASTER CURVE	
		g/L	sec
Optical	Yumizen G400 DDi/G200	4	9,0
		3	11,0
		2	14,5
		1	23,5
		Dilution	10x



06/05/2021