



# Yumizen G FIB 2 for Yumizen G1500/G1550 and Yumizen G800

<b>Reference No.:</b>		<b>1300036383</b>	
<b>Lot No.:</b>		<b>921110 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2024-11-30</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
		<b>g/L</b>	<b>sec</b>
<b>Optical</b>	Yumizen G1500/G1550/G800	4	8,6
		3	10,3
		2	13,4
		1	21,0
		<b>Dilution</b>	<b>10x</b>



# Yumizen G FIB 2 for Yumizen G400 DDi/G200

<b>Reference No.:</b>		<b>1300036383</b>	
<b>Lot No.:</b>		<b>921110 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2024-11-30</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
		<b>g/L</b>	<b>sec</b>
<b>Optical</b>	Yumizen G400 DDi/G200	4	9,4
		3	11,5
		2	15,3
		1	24,9
		<b>Dilution</b>	<b>10x</b>



06/05/2021